

Re: Olympus BH-2

Source: <http://sci.tech-archive.net/Archive/sci.techniques.microscopy/2005-04/msg00025.html>

- *From:* "GTO" <gregor_o@xxxxxxxxxxxxxxxxxxxx>
 - *Date:* Tue, 05 Apr 2005 22:10:32 GMT
-

Trust me, after spending 4 to 6k on a cooled CCD imaging device, spending another 1k for a decent mount that is compatible is most likely not an issue.

Check out <http://www.diaginc.com/> and their adaptaters. One of their imaging devices is on my wishlist ;-)

Gregor

"Andrew" <andrewmoffat@xxxxxxxx> wrote in message
news:1112703350.116654.151940@xx
>I have on Olympus BH-2 microscope fitted with a 1.4 MegaPixel DVC
> digital camera. I use the microscope for materials analysis (mostly on
> metals and some polymers). I am thinking of upgrading my digital camera
> to a colour one but I am not sure of the best one to get so I am after
> some advice please.
>
> One of the most important things for the new camera would be that it
> has good quality square pixels so that I can measure the grains etc...
> accurately, I don't know if this is a problem with most newer cameras.
> The resolution is fine at 1.4 megapixels. I do use polarised light so I
> think that peltier cooled is better and I would preferably like not to
> have to pay for new adaptors so a standard C mount fit with a 2/3 inch
> ccd would be good.
>
> Any advice will be gratefully received. Thanks,
>
> Andrew
>

-
- *References:*
 - ◆ *Olympus BH-2*
 - ◇ *From:* Andrew

Re: Olympus BH-2

- Prev by Date: ***Re: Magnification factor of my digital camera***
- Next by Date: ***Turbo Pump and Ion Pump Repairs.***
- Previous by thread: ***Olympus BH-2***
- Next by thread: ***ES: EG&G Photometric Microscope System***
- Index(es):
 - ◆ ***Date***
 - ◆ ***Thread***