

## Re: High NA darkfield condenser

---

*Source:* <http://sci.tech-archive.net/Archive/sci.techniques.microscopy/2006-02/msg00049.html>

---

- *From:* "GTO" <[gregor\\_o@xxxxxxxxxxxxxxxx](mailto:gregor_o@xxxxxxxxxxxxxxxx)>
  - *Date:* Mon, 13 Feb 2006 08:23:22 GMT
- 

Good point. This should work for some types of mounts (eg. permanent mounts).

Gregor

"Aaron" <[ngly@xxxxxxxx](mailto:ngly@xxxxxxxx)> wrote in message  
[news:pm2ru1lcdofsl2eigaice97trfq4rj661f@xxxxxxxx](mailto:news:pm2ru1lcdofsl2eigaice97trfq4rj661f@xxxxxxxx)

Hi Gregor,

The answer is yes. You just invert the slide and view it coverglass down.

Aaron

On 10 Feb 2006 11:43:35 -0800, "GTO" <[Gregor.Overney@xxxxxxxx](mailto:Gregor.Overney@xxxxxxxx)> wrote:

Can you oil a condenser to a dish/slide with an inverted scope?

Gregor