

## Re: How to setup Zeiss DIC?

---

*Source:* <http://sci.tech--archive.net/Archive/sci.techniques.microscopy/2008-07/msg00015.html>

---

- *From:* Perl <thl438v@xxxxxxxxxx>
  - *Date:* Mon, 7 Jul 2008 07:24:59 -0700 (PDT)
- 

Hi

You need following:

Polarizer  
INKO Condenser (with DIC Filter)  
Plan or better Objectives  
DIC Slider matching Objectives in Holder (Zeiss)  
Analyser for DIC marked II or III in you case

Alignment and set up instruction can i mail you  
if you send you email

PS / There will take some time to setup first time and  
a Phase telescope will help a lot

Best Regards  
Pär Lundqvist  
Sweden

.