

Re: adaptaters for Leica DM 500

Source: <http://sci.tech--archive.net/Archive/sci.techniques.microscopy/2009-06/msg00063.html>

- *From:* David <pennine56@xxxxxxxxxxxx>
 - *Date:* Mon, 15 Jun 2009 09:20:07 -0700 (PDT)
-

On Jun 15, 11:24 am, Daniel <danieldezut...@xxxxxxxxxx> wrote:

Hello,
I bought a Leica DM 500 microscope with a trinocular head (inlet diameter: 37 mm).
I would like to know what adapters do I need to place a mini digital camera on the outlet of the trinocular microscope?
What other adapters do I need to place a camera canon eos 400 D?
Thank you

Daniel

Hello Daniel

The manual with systems diagram is available to download from various sites if you do a Google search. The manual shows that the trino head is a C mount for appropriate relay lens so depending on budget you can explore options that suit either Leica's relay lenses and cameras or elsewhere. Cheapest may be a security style camera with a video grabber or could choose one of the many microscope cameras with C mount.

There are C mount to T2 SLR camera lens adapters if the the microscope fitting is female but may not be a good match of image size with DSLR sensor size, you could play around with tracing paper to inspect projected image to get a feel for what may work.

David

.